
UpFuse - Characterization

28H4

Jul 10, 2020

CONTENTS:

| | | |
|----------|------------------------|----------|
| 1 | Requirements | 1 |
| 2 | measurement.py | 2 |
| 3 | reset_smu.py | 3 |
| 4 | Keithley236 API | 4 |
| | Index | 5 |

REQUIREMENTS

- The script is based on PyVisa. This package requires a VISA interface in the background ([NI-VISA download via: ni.com](#))
- Driver for the GPIB-USB-interface ([NI-488.2 download via: ni.com](#))

MEASUREMENT.PY

Script to perform measurements for the characterization of nanoionic devices.

`measurement.measurement` (*file_set_times*, ***kwargs*)

Performs measurements for the characterization of nanoionic memories using a Keihtley 236.

Parameters

:param file_set_times: str The file location of an Excel file containing the individual set times.

:key gpib_address: int GPIB adress of the Keihtley 236. (Default is 16).

:key compliance: str Compliance Value [A]. Scientific notation required e.g “1E-9” for 1nA. (Default value is 1E-9).

:key measurement_delay: float Delay time [s] between applying the measuring voltage and the start of the measurement. (Default value is 30).

:key measurement_range: str Measurement source_range for the current. Allowed values are: Auto, 1nA, 10nA, 100nA, 1μA, 10μA, 100μA, 1mA, 10mA and 100mA. (Default value is 1nA).

:key measurement_voltage: float Applied voltage [V] during the measurement step. (Default value is 1).

:key rest_period: float Delay time [s] between the end of the measurement and the start of the next set pulse. (Default value is 120).

:key set_voltage: float Applied voltage [V] during the set step. (Default value is 1).

`measurement.read_set_times` (*file*)

Returns the first column of an Excel file as list.

RESET_SMU.PY

Simple script to reset the Keithley236 to its factory defaults.

KEITHLEY236 API

```
class SimpleKeithley236.Keithley236.Keithley236(gpib_address=16,  
                                              compliance_level='1E-3',    measure-  
                                              ment_range='1mA')  
    Minimal control interface (GPIB) for a Keithley 236 for the characterization of nano-ionic memories.
```

Methods

| | |
|--|---|
| <code>impulse(voltage, duration)</code> | Applies the voltage [V] for the duration [s]. |
| <code>measurement(voltage[, delay])</code> | Returns measured current [A] for the applied voltage [V]. |

INDEX

K

Keithley236 (class in SimpleKeith-
ley236.Keithley236), 4

M

measurement
module, 2

measurement() (in module measurement), 2

module
measurement, 2
reset_smu, 3

R

read_set_times() (in module measurement), 2

reset_smu
module, 3